

Lymphozyten (rel.)

Einheit	%																		
Häufigkeit	werktags, normalschichtig																		
Material	EDTA-Blut																		
Probengefäß	EDTA-Monovette																		
Probenvolumen	2.7 ml																		
Qualitätsmanagement-Status	akkreditiert																		
Referenzbereich	<table border="1"><thead><tr><th>Alter</th><th>Referenzbereich</th></tr></thead><tbody><tr><td>≤ 5T</td><td>22,0 - 37,0</td></tr><tr><td>5T - 1W6T</td><td>40,0 - 81,0</td></tr><tr><td>1W6T - 1J</td><td>40,0 - 85,0</td></tr><tr><td>1J - 3J</td><td>50,0 - 70,0</td></tr><tr><td>3J - 5J</td><td>37,0 - 52,0</td></tr><tr><td>5J - 9J</td><td>30,0 - 48,0</td></tr><tr><td>9J - 16J</td><td>33,0 - 48,0</td></tr><tr><td>> 16J</td><td>20,0 - 49,0</td></tr></tbody></table>	Alter	Referenzbereich	≤ 5T	22,0 - 37,0	5T - 1W6T	40,0 - 81,0	1W6T - 1J	40,0 - 85,0	1J - 3J	50,0 - 70,0	3J - 5J	37,0 - 52,0	5J - 9J	30,0 - 48,0	9J - 16J	33,0 - 48,0	> 16J	20,0 - 49,0
Alter	Referenzbereich																		
≤ 5T	22,0 - 37,0																		
5T - 1W6T	40,0 - 81,0																		
1W6T - 1J	40,0 - 85,0																		
1J - 3J	50,0 - 70,0																		
3J - 5J	37,0 - 52,0																		
5J - 9J	30,0 - 48,0																		
9J - 16J	33,0 - 48,0																		
> 16J	20,0 - 49,0																		
Verfahren	Elektronische Zählung																		
letzte Änderung	29.11.2024 11:57																		