

Neutr. Granulozyten (rel.)

Abnahmevorschrift																			
Einheit	%																		
Häufigkeit	werktags, normalschichtig																		
Material	EDTA-Blut																		
Probengefäß	EDTA-Monovette																		
Probenvolumen	2.7 ml																		
Qualitätsmanagement-Status	akkreditiert																		
Referenzbereich	<table><thead><tr><th>Alter</th><th>Referenzbereich</th></tr></thead><tbody><tr><td>≤ 5T</td><td>67,0 - 90,0</td></tr><tr><td>5T - 1W6T</td><td>31,0 - 48,0</td></tr><tr><td>1W6T - 1J</td><td>38,0 - 49,0</td></tr><tr><td>1J - 3J</td><td>26,0 - 50,0</td></tr><tr><td>3J - 5J</td><td>30,0 - 53,0</td></tr><tr><td>5J - 9J</td><td>33,0 - 53,0</td></tr><tr><td>9J - 16J</td><td>43,0 - 57,0</td></tr><tr><td>> 16J</td><td>36,0 - 77,0</td></tr></tbody></table>	Alter	Referenzbereich	≤ 5T	67,0 - 90,0	5T - 1W6T	31,0 - 48,0	1W6T - 1J	38,0 - 49,0	1J - 3J	26,0 - 50,0	3J - 5J	30,0 - 53,0	5J - 9J	33,0 - 53,0	9J - 16J	43,0 - 57,0	> 16J	36,0 - 77,0
Alter	Referenzbereich																		
≤ 5T	67,0 - 90,0																		
5T - 1W6T	31,0 - 48,0																		
1W6T - 1J	38,0 - 49,0																		
1J - 3J	26,0 - 50,0																		
3J - 5J	30,0 - 53,0																		
5J - 9J	33,0 - 53,0																		
9J - 16J	43,0 - 57,0																		
> 16J	36,0 - 77,0																		
Verfahren	Elektronische Zählung																		
letzte Änderung	02.01.2025 00:31																		