

Eosino. Granulozyten (rel.)

Einheit	%																	
Häufigkeit	werktags, normalschichtig																	
Material	EDTA-Blut																	
Probengefäß	EDTA-Monovette																	
Probenvolumen	2.7 ml																	
Qualitätsmanagement-Status	akkreditiert																	
Referenzbereich	<table><thead><tr><th>Alter</th><th>Referenzbereich</th></tr></thead><tbody><tr><td>≤ 5T</td><td>0,0 - 3,0</td></tr><tr><td>5T - 1J</td><td>1,0 - 5,0</td></tr><tr><td>1J - 3J</td><td>1,0 - 4,0</td></tr><tr><td>3J - 5J</td><td>0,0 - 4,0</td></tr><tr><td>5J - 9J</td><td>0,0 - 5,0</td></tr><tr><td>9J - 16J</td><td>0,0 - 4,0</td></tr><tr><td>> 16J</td><td>0,0 - 5,0</td></tr></tbody></table>		Alter	Referenzbereich	≤ 5T	0,0 - 3,0	5T - 1J	1,0 - 5,0	1J - 3J	1,0 - 4,0	3J - 5J	0,0 - 4,0	5J - 9J	0,0 - 5,0	9J - 16J	0,0 - 4,0	> 16J	0,0 - 5,0
Alter	Referenzbereich																	
≤ 5T	0,0 - 3,0																	
5T - 1J	1,0 - 5,0																	
1J - 3J	1,0 - 4,0																	
3J - 5J	0,0 - 4,0																	
5J - 9J	0,0 - 5,0																	
9J - 16J	0,0 - 4,0																	
> 16J	0,0 - 5,0																	
Verfahren	Elektronische Zählung																	
letzte Änderung	29.11.2024 11:30																	