

Monozyten (rel.) im EDTA-Blut (%)

Einheit	%																		
Häufigkeit	werktags, normalschichtig																		
Material	EDTA-Blut																		
Probengefäß	EDTA-Monovette																		
Probenvolumen	2.7 ml																		
Qualitätsmanagement-Status	akkreditiert																		
Referenzbereich	<table border="1"><thead><tr><th>Alter</th><th>Referenzbereich</th></tr></thead><tbody><tr><td>≤ 5T</td><td>4,0 - 10,0</td></tr><tr><td>5T - 1W6T</td><td>6,0 - 13,0</td></tr><tr><td>1W6T - 1J</td><td>4,0 - 12,0</td></tr><tr><td>1J - 3J</td><td>1,0 - 6,0</td></tr><tr><td>3J - 3J</td><td>1,0 - 6,0</td></tr><tr><td>3J - 5J</td><td>0,0 - 5,0</td></tr><tr><td>5J - 16J</td><td>0,0 - 6,0</td></tr><tr><td>> 16J</td><td>0,0 - 9,0</td></tr></tbody></table>	Alter	Referenzbereich	≤ 5T	4,0 - 10,0	5T - 1W6T	6,0 - 13,0	1W6T - 1J	4,0 - 12,0	1J - 3J	1,0 - 6,0	3J - 3J	1,0 - 6,0	3J - 5J	0,0 - 5,0	5J - 16J	0,0 - 6,0	> 16J	0,0 - 9,0
Alter	Referenzbereich																		
≤ 5T	4,0 - 10,0																		
5T - 1W6T	6,0 - 13,0																		
1W6T - 1J	4,0 - 12,0																		
1J - 3J	1,0 - 6,0																		
3J - 3J	1,0 - 6,0																		
3J - 5J	0,0 - 5,0																		
5J - 16J	0,0 - 6,0																		
> 16J	0,0 - 9,0																		
Verfahren	mikroskopische Zählung																		
letzte Änderung	29.11.2024 11:57																		